

NeuroTherm Part Number	Electrode Compatibility	Generator Compatibility
DAC-NT	NeuroTherm Disposable Electrode	NeuroTherm Generator
DACUK-NT	NeuroTherm International Disposable Electrode	NeuroTherm Generator
AC-NT-C	Reusable Cordotomy Electrode	NeuroTherm Generator
AC-NT-COS	Cosman Reusable Electrode	NeuroTherm Generator
AC-NT-SN	Smith & Nephew Reusable Electrode	NeuroTherm Generator
AC-NT-TC2	Radionics TC2 Electrode	NeuroTherm Generator
AC-Bi-P	NeuroTherm Reusable Electrode	NeuroTherm Generator
AC-TC2-NT	NeuroTherm Reusable Electrode	Radionics Generator
DAC-RAD	NeuroTherm Disposable Electrode	Radionics Generator
DACUK-RAD	NeuroTherm International Disposable Electrode	Radionics Generator
AC-COS-NT	NeuroTherm Reusable Electrode	Cosman Generator
DAC-COS	NeuroTherm Disposable Electrode	Cosman Generator
DACUK-COS	NeuroTherm International Disposable Electrode	Cosman Generator
AC-B-NT	NeuroTherm Reusable Electrode	Kimberly Clark/Baylis Generator
AC-B-DUS	NeuroTherm Disposable Electrode	Kimberly Clark/Baylis Generator
AC-B-DUK	NeuroTherm International Disposable Electrode	Kimberly Clark/Baylis Generator
AC-SN-NT	NeuroTherm Reusable Electrode	Smith & Nephew Generator
AC-SN-DUS	NeuroTherm Disposable Electrode	Smith & Nephew Generator
AC-SN-DUK	NeuroTherm International Disposable Electrode	Smith & Nephew Generator
DACUK-NT-L	NeuroTherm Disposable Electrode	NeuroTherm Generator
AC-SI	NeuroTherm Disposable Simplicity Electrode	NeuroTherm Generator
AC-IDET-4	Smith & Nephew Disposable IDET Cath (4 Pin)	NeuroTherm Generator
AC-CN-S	SMK Radionics Electrode	NeuroTherm Generator
AC-CN-R	RFK Radionics Electrode	NeuroTherm Generator
AC-CN-C	Radionics Cordotomy Electrode	NeuroTherm Generator
AC-POLE-NT	Disposable Pole Needle	NeuroTherm Generator
AC-C-S	SMK Radionics Electrode	Radionics Generator
AC-C-R	RFK Radionics Electrode	Radionics Generator
AC-POLE-RAD	Disposable Pole Needle	Radionics Generator